


Non-Woven Single-Coated Tape



Backing:

Adhesive:

1.3 oz/sy white spunlace polyester nonwoven

Non-sensitizing medical-grade acrylic (M109)

Liner:

72lb white polycoated kraft, silicone treated on one side

3614B

- Electrodes: ECG, EKG, TENS
- Wound coverings

- Applications
 - Medical device attachment
 - Island dressings

Features	Benefits
PET non-woven carrier	Breathable for maximum patient comfort
	Highly conformable
 Medical-grade Acrylic Adhesive 	 Designed for short-term and medium wear uses
	 Stable to EtO and gamma sterilization
72lb Release Liner	Easy to die-cut and convert
	 Good dimensional stability and lay flat

Typical Test Values

Metric UOM	Test Method
0.218 mm	PSTC 133
0.051 mm	PSTC 133
0.142 mm	PSTC 133
8.75 N/cm	PSTC 101 Mod
	0.218 mm 0.051 mm 0.142 mm

The base adhesives in this medical product have been tested in accordance with ISO 10993 and found to be noncytotoxic, non-irritating, and non-sensitizing.

In accordance with FDA cGMP requirements, Berry reminds the customer that ultimate fitness for use must be determined by the manufacturer of the finished device.

Recommended Storage Conditions: 40-60% Humidity, 60°F-80°F

The information presented herein was prepared at BERRY GLOBAL, INC. (25 Forge Parkway Franklin, MA 02038) by qualified technical personnel. To BERRY's knowledge it is true and accurate. However, the information and recommendations are furnished for these products with the understanding that the purchaser or user, as the case may be, will independently determine that the product is suitable for the intended use, and that such use complies with all applicable federal, state and local laws and regulations. The data are submitted only for the user's information and consideration, and do not constitute a warranty of any kind (including but not limited to a warranty of merchantability or of fitness for a particular purpose), or a representation for which BERRY assumes any legal responsibility.



25 Forge Pkwy Franklin, MA 02038 1-800-343-7875 Fax 1-800-328-4822 www.berryglobal.com